

 **SERVImed**
 **iROMED Group**

CXL

CORNEAL CROSS-LINKING

Y
CARO



SERVIMED'S NEW AND IMPROVED CROSS-LINKING SYSTEM

VERSATILE OPEN SYSTEM

- 10" TOUCHSCREEN DISPLAY
- INTEGRATED CAMERA
- UV POWER: 1-30 mW/cm²
- CONTINUOUS AND PULSED
- 3-12 mm ADJUSTABLE DIAMETER
- LOCKING WHEELS
- MULTI-POINT ARM LOCKING

WIDE TREATMENT RANGE

- DRESDEN EPI-OFF; STANDARD EPI-ON;
- IONTOPHORESIS; ACCELERATED; USER
- PRESETS; KERATITIS; CUSTOM FAST CXL*
- (*only with RIBOCROSS te[®] and RIBOFAST)

ALWAYS CONNECTED

- LAN AND WIFI READY
- TREATMENT REPORT ARCHIVE
- USB STORAGE AND PRINT
- REMOTE ASSISTANCE

OUR EXCLUSIVE CUSTOMIZED CXL TREATMENT

CUSTOMIZED

- Topo-pachimetric guided treatment
- Customized procedure parameters
- Non invasive EPI-ON technique
- Endothelial safeguard with no additional soaking needed
- Suitable for thin corneas

INNOVATIVE

- Vitamin E TPGS – enhanced riboflavin for improved corneal penetration with safe and effective results

PEER-REVIEWED

- Multiple peer-reviewed publications [1-8]
- 7-year follow-up published on Cornea [7]





SERVImed Industrial S.p.A.

Via Tempio del Cielo 3/5, 00144 Rome (Italy)
Tel: +39 06 92595490 Fax: +39 06 89360010
Email: info@servimedindustrial.com
www.servimed-industrial.com



Iromed Group S.r.l.

Via Tempio del Cielo 3/5, 00144 Rome (Italy)
Email: info@iromedgroup.com
www.iromedgroup.com